

TORRANCE MEMORIAL MICROBIOLOGY SWAB GUIDE- OUTREACH

- **Anaerobic/Aerobic Culture**
 - Test Code 63060235
- **Aerobic Culture**
 - Wounds, skin, etc.
 - Test Code 63060323
- **Group A Strep PCR Throat**
 - Test Code 2326274037
- **MRSA Evaluation**
 - Test Code 95259289
- **Throat Culture**
 - Test Code 63060297
- **Group B Strep Culture**
 - Test Code 95259283

Bacterial Swab

Use ESwab



- **SARS-CoV-2/Flu/RSV PCR**
 - Test Code 2080288029
- **Respiratory PCR Panel**
 - Test Code 394064397
- **SARS-CoV-2 PCR**
 - NP, Nasal, Throat
 - Test Code 1785398229
- **Herpes Simplex Virus 1-2/Varicella Zoster Virus PCR (non-Blood or CSF)**
 - Test Code 6408427429

Viral Swab

Use Universal (Viral) Transport Media (UTM or UVT)
Pediatric (1mL) or Adult (3mL)



- **Chlamydia trachomatis/Neisseria gonorrhoeae/Trichomonas vaginalis PCR**
 - Test Code 4009603791
- **Vaginal PCR Panel**
 - Test Code 947198625

**Transfer urine from cup
IMMEDIATELY.**

SEND TWO SWABS if both
CT/GC/TV and Vaginal PCR panels
are ordered.

BD Molecular Buffer Tube

**Urine
Vaginal**



- **Chlamydia trachomatis and Neisseria gonorrhoeae, Nucleic Acid Amplification, Varies**
 - Test Code: 256128714

ORANGE LABEL SOURCES:

Throat/Oropharynx
Rectal/Anal
Vaginal (if no **green** BD)

WHITE LABEL SOURCES:

Urethra (Males)
Cervix/Endocervix (if no **green** BD)

Hologic Aptima Swab

Use Hologic Aptima
**Oral
Rectal**
Non-Urine or Non-Female Genital Sources

